

## Release Note for Sense notetaker version 7.0 software patch

### New Features

#### 1. Sense Chat

Available in “Social Networking” menu, used for chatting between Sense notetakers. Useful for any 2 people with Sense notetakers, but particularly useful for Deaf-Blind people to use 2 Sense notetakers to talk to each other.

- **Connection**  
When you open Sense chat, choose “Bluetooth Serial Send” on one device and “Bluetooth Serial Receive” on the other device. One will ask for a pin code, then the other. Enter the same pin code on both devices. Conversation will then be available between the 2 devices.

See section 16.5 of the User Manual for more details. Looks very similar to Google Talk dialog as outlined below:

- **Main Control**  
“Send Message” edit box, “Receive Message” edit box, Close button
- **Option settings:**
  - a) Message Alert Setting: Alert and Display/Alert only
  - b) Send/Receive Message Alert: sound1/sound2/vibration
  - c) Alert when messages are sent/received/when users connect/disconnect: Sound/Vibration/Off

#### 2. Google Search

- **Controls:** “Search Word” edit box, “Search” button, “Close” button, “Category” combo box
- Located in “Extras” menu, see section 15.4 of the User Manual

#### 3. Book Share Download

- Institutional accounts are now supported

#### 4. Web Browser

- Move to previous visited link: Space-dots-2-3-6 (Alt-Shift-v)
- Move to next visited link: Space-dots-3-5-6 (Alt-v)

#### 5. Backup and Restore

- Inserted back-up of web browser cookies

## 6. Global Options

- Braille Code(c)
  - a) Added the ability to use U.K Braille

## Changes/Features

### 1. General

- Sense Navigation is moved to “Extra’s” menu
- Changed the keyboard layout for QWERTY products to English International
- Now remembers Braille grade when switching Braille codes
- Added hot-key to change the “Sub Voice Rate”: (Braille: Enter-Space-F1 / Enter-Space-F4, Qwerty: Shift-Windows-F6/ Windows-F6)

### 2. Word Processor

- Inserted the option to select whether you want the notetaker to load formatting information (Off by default)

Access in Settings → Support formatting documents (f) option

- Inserted the “Read Selected Text in Current Language” option in the “Read” menu
- Changed the QWERTY hot-key to input control characters from FN-Space to Ctrl-Space
- Changed other hot-keys which formerly used Ctrl-Alt
  - a) Select edit language attribute: Alt-Windows-i
  - b) Select change language: Alt-Windows-v
  - c) Read selected text with current language: Alt-Windows-r

### 3. Web Browser

- Inserted the “Go to Previous Visited Link” and “Go to Next Visited Link” items in the “Go To” menu
- Now ignores security alert dialog boxes on all websites

### 4. Changed the Spanish Translation Rules

- á (a acute): Dots 1-2-3-4-5
- é (e acute): Dots-2-3-4-6
- í (i acute): Dots-3-4
- ñ (n tilde): Space-U, Dots-1-2-4-5-6
- ó (o acute): Dots-3-4-6

- ú (u acute): Dots-2-3-4-5-6
- ç (c cedilla): Dots-1-2-3-4-6

## 5. Sense Dictionary

- Changed QWERTY hot-keys.
  - a) English Dictionary: Ctrl-Windows-d
  - b) English Thesaurus: Ctrl-Windows-t
  - c) Spanish Dictionary: Ctrl-Windows-l
  - d) French Dictionary: Ctrl-Windows-f
  - e) Italian Dictionary: Ctrl-Windows-y
- Improved the function to search words ending in “ies”

## 6. Utilities Menu Calendar

- No longer changes calendars when pressing Enter after typing date.

## 7. Options

- Changed the default setting for “Punctuation Level” from “all” to “off”
- Changed the “Braille Code” hot-key from Ctrl-Alt-b to Alt-Windows-b

## Bug Fixes

### 1. General

- Fixed the problem of incorrect cursor routing

### 2. Word Processor

- Fixed the problem of accessing table information when formatting is set to on
- Fixed the problem of crashing when opening special doc documents
- Fixed the problem of printing format characters when embossing or pasting to another program such as e-mail
- Fixed the issue of deleting data when saving documents
- Fixed sluggishness when opening/saving Word documents

### 3. Daisy Player

- Fixed the problem of not displaying all data in Braille

### 4. Terminal for Screen Reader

- Fixed the problem of inability to connect with Mountain Lion using Bluetooth (Braille Sense PLUS and QWERTY)
- Fixed the problem of not announcing the connection options when the “Announce Shortcuts” option was turned off

#### **5. Schedule Manager**

- Fixed problems restoring the schedule data (Braille Sense U2 and QWERTY)

#### **6. E-mail**

- Fixed the problem of not activating the “Modify” and “Cancel” button in the “Set Path” dialog box
- Fixed the problem of mistranslated e-mail headers on reply and forward messages
- Fixed the problem of displaying incorrect data
- Fixed the problem using the scroll keys in the advance dialog box of the Accounts Manager
- Fixed inability to create new mail box on first open
- Fixed inability to navigate the mail box list using first letter navigation
- Fixed problems opening HTML email
- Fixed the problem of displaying incorrect prompts

#### **7. Radio**

- Fixed the problem of not displaying the frequency after closing the “Settings” dialog box

#### **8. Global Options**

- Fixed inability to use the scroll keys when using “typing mode”

#### **9. Media player**

- Fixed the problem of crashing after playing just a few songs when using shuffle mode

#### **10. Backup/Restore Personalized Settings**

- Fixed the problem of not restoring settings correctly
- Now remembers switch Control/Function and Alt/Command for U2 QWERTY

#### **11. Calendar**

- Fixed the problem of not being able to move between Lunar and Solar calendars using the tab and shift-tab keys